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The screenshot shows a software interface with a menu bar at the top: File, Reports, Games, Setup, Language, Help. The 'Reports' menu is open, showing options: Progress Report, Home Report, Device Report, and Printouts and Export. The 'Printouts and Export' sub-menu is also open, showing: Export Patient List To CSV, Export all session statistics to CSV, Print patient details, Print patient sessions, Print patient list, Print patient session statistics, and Open Exported Data. The 'Open Exported Data' option is highlighted. Below the menu, there is a table with columns: PatID, Surname, First name, Protocol, Condition, Attachments, and Notes. The table contains 10 rows of patient session data. Below the table, there are checkboxes for 'Confidential' and 'Default Patient'. At the bottom, there are two sections: 'Patient Options' and 'Session Options', each with several buttons for actions like 'View / Edit', 'View Device Report', 'Add', 'Delete', 'Review', 'Delete Session', 'Deselect all sessions', 'Create Progress report (Assessment)', 'Start Protocol', and 'Backup data'.

PatID	Surname	First name	Protocol	Condition	Attachments	Notes
0~	Example	Bob				
1	Nunes	Anna	Device Controlled ETS	Multiple Sclerosis	Attachments (0)	
2	Nuggets	Robert	Open Display	Neurofeedback	Attachments (0)	spasticity
3	Menjega	Alex	Open Display	Neurofeedback	Attachments (0)	after stim
4	Stewart	Alice	Device Controlled EMG	Neurofeedback	Attachments (0)	
5	Moss	Maurice	Device Controlled ETS	Neurofeedback	Attachments (0)	
6			Open Display	Neurofeedback	Attachments (0)	
7			Device Controlled STIM	Neurofeedback	Attachments (0)	
8			Device Controlled EMG	Neurofeedback	Attachments (0)	
9			Device Controlled EMG	Neurofeedback	Attachments (0)	
10			Open Display	Neurofeedback	Attachments (0)	

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